

Product Datasheet

Mouse Protein-glutamine Gamma-glutamyltransferase K ELISA Kit (orb1210404)

Description	Mouse Protein-glutamine Gamma-glutamyltransferase K, TGM1 ELISA Kit
Reactivity	Mouse
Range	0.5-300mU/mL
Tested Applications	ELISA
Target	TGM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	89,826
Uniprot ID	Q9JLF6
Solubility (25°C)	320mU/mL
Sensitivity	0.27mU/mL
Entrez	21816
Expiration Date	Please enquire.