

Product Datasheet

Mouse Protein-arginine Deiminase Type-2 ELISA Kit (orb1210405)

Description	Mouse Protein-arginine Deiminase Type-2, PADI2 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	PADI2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	76,250
Uniprot ID	Q08642
Solubility (25°C)	6400ng/L
Sensitivity	9.55ng/L
Entrez	18600
Expiration Date	Please enquire.