

Product Datasheet

Mouse Protein Z-dependent Protease Inhibitor ELISA Kit (orb1210406)

Description	Mouse Protein Z-dependent Protease Inhibitor, SERPINA10 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	SERPINA10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,797
Uniprot ID	Q8R121
Solubility (25°C)	4800ng/L
Sensitivity	99.61ng/L
Entrez	217847
Expiration Date	Please enquire.