

Product Datasheet

Mouse Protein Wnt-8B ELISA Kit (orb1210407)

Description	Mouse Protein Wnt-8B, WNT8B ELISA Kit
Reactivity	Mouse
Range	20-3800ng/L
Tested Applications	ELISA
Target	WNT8B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,576
Uniprot ID	Q9WUD6
Solubility (25°C)	8000ng/L
Sensitivity	8.85ng/L
Entrez	22423
Expiration Date	Please enquire.