

Product Datasheet

Mouse Protein Wnt-8A ELISA Kit (orb1210408)

Description	Mouse Protein Wnt-8A, WNT8A ELISA Kit
Reactivity	Mouse
Range	20-7000ng/L
Tested Applications	ELISA
Target	WNT8A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,573
Uniprot ID	Q64527
Solubility (25°C)	9600ng/L
Sensitivity	10.52ng/L
Entrez	20890
Expiration Date	Please enquire.