

Product Datasheet

Mouse Protein Wnt-16 ELISA Kit (orb1210413)

Description	Mouse Protein Wnt-16, WNT16 ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	WNT16
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,727
Uniprot ID	Q9QYS1
Solubility (25°C)	9600ng/L
Sensitivity	14.06ng/L
Entrez	93735
Expiration Date	Please enquire.