

Product Datasheet

Mouse Protein Notum Homolog ELISA Kit (orb1210418)

Description	Mouse Protein Notum Homolog, NOTUM ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	NOTUM
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,005
Uniprot ID	Q0VBF8
Solubility (25°C)	24ng/mL
Sensitivity	0.019ng/mL
Entrez	381310
Expiration Date	Please enquire.