

Product Datasheet

Human Protein Wnt-5A ELISA Kit (orb1210448)

Description	Human Protein Wnt-5A, WNT5A ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	WNT5A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,339
Uniprot ID	P41221
Solubility (25°C)	32ng/ml
Sensitivity	0.027ng/ml
Entrez	7474
Expiration Date	Please enquire.