

Product Datasheet

Human Protein Wnt-4 ELISA Kit (orb1210449)

Description	Human Protein Wnt-4, WNT4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	WNT4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,052
Uniprot ID	P56705
Solubility (25°C)	3200ng/L
Sensitivity	7.66ng/L
Entrez	54361
Expiration Date	Please enquire.