

Product Datasheet

Mouse Phenylalanine-4-hydroxylase ELISA Kit (orb1210457)

Description	Mouse Phenylalanine-4-hydroxylase, PAH ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	PAH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,900
Uniprot ID	P16331
Solubility (25°C)	3200ng/L
Sensitivity	7.92ng/L
Entrez	18478
Expiration Date	Please enquire.