

Product Datasheet

Human Usp6 N-terminal Like Protein ELISA Kit (orb1210458)

Description	Human Usp6 N-terminal Like Protein, USP6NL ELISA Kit
Reactivity	Human
Range	20-3500ng/L
Tested Applications	ELISA
Target	USP6NL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	94,104
Uniprot ID	Q92738
Solubility (25°C)	4000ng/L
Sensitivity	9.35ng/L
Entrez	9712
Expiration Date	Please enquire.