

Product Datasheet

Rat VGF Nerve Growth Factor Inducible ELISA Kit (orb1210478)

Description	Rat VGF Nerve Growth Factor Inducible, VGF ELISA Kit
Reactivity	Rat
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	VGF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,010
Uniprot ID	P25427
Solubility (25°C)	48ng/mL
Sensitivity	0.053ng/mL
Entrez	310738
Expiration Date	Please enquire.