

Product Datasheet

Human Claudin 11 ELISA Kit (orb1210502)

Description	Human Claudin 11, CLDN11 ELISA Kit
Reactivity	Human
Range	1-40ng/ml
Tested Applications	ELISA
Target	CLDN11
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,993
Uniprot ID	075508
Solubility (25°C)	48ng/ml
Sensitivity	0.46ng/ml
Entrez	5010
Expiration Date	Please enquire.