

Product Datasheet

Rat Troponin I ELISA Kit (orb1210508)

Description	Rat Troponin I, Slow Skeletal Muscle, TNNI1 ELISA Kit
Reactivity	Rat
Range	6-1800ng/L
Tested Applications	ELISA
Target	TNNI1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,724
Uniprot ID	P13413
Solubility (25°C)	1920ng/L
Sensitivity	3.42ng/L
Entrez	29388
Expiration Date	Please enquire.