

Product Datasheet

Human Transgelin 2 ELISA Kit (orb1210527)

Description	Human Transgelin 2, TAGLN2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	TAGLN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,391
Uniprot ID	P37802
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	8407
Expiration Date	Please enquire.