

Product Datasheet

Mouse Muscarinic Acetylcholine Receptor M2 ELISA Kit (orb1210549)

Description	Mouse Muscarinic Acetylcholine Receptor M2, CHRM2 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	CHRM2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,499
Uniprot ID	Q9ERZ4
Solubility (25°C)	48ng/mL
Sensitivity	0.059ng/mL
Entrez	243764
Expiration Date	Please enquire.