

Product Datasheet

Rat Serine Threonine-protein Kinase Tnni3K ELISA Kit (orb1210574)

Description	Rat Serine Threonine-protein Kinase Tnni3K, TNNI3K ELISA Kit
Reactivity	Rat
Range	10-2000ng/L
Tested Applications	ELISA
Target	TNNI3K
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	92,732
Uniprot ID	Q7TQP6
Solubility (25°C)	2400ng/L
Sensitivity	4.79ng/L
Entrez	295531
Expiration Date	Please enquire.