

Product Datasheet

Rat Runt-related Transcription Factor 2 ELISA Kit (orb1210586)

Description	Rat Runt-related Transcription Factor 2, RUNX2 ELISA Kit
Reactivity	Rat
Range	1-400ng/mL
Tested Applications	ELISA
Target	RUNX2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,247
Uniprot ID	Q9Z2J9
Solubility (25°C)	480ng/mL
Sensitivity	0.46ng/mL
Entrez	367218
Expiration Date	Please enquire.