

Product Datasheet

Human Bone Morphogenetic Protein 5 ELISA Kit (orb1210587)

Description	Human Bone Morphogenetic Protein 5, BMP5 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	BMP5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,737
Uniprot ID	P22003
Solubility (25°C)	2400ng/L
Sensitivity	4.88ng/L
Entrez	653
Expiration Date	Please enquire.