

Product Datasheet

Mouse Matrix Metalloproteinase-14 ELISA Kit (orb1210598)

Description	Mouse Matrix Metalloproteinase-14, MMP14 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	MMP14
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	65,919
Uniprot ID	P53690
Solubility (25°C)	24ng/mL
Sensitivity	0.022ng/mL
Entrez	17387
Expiration Date	Please enquire.