

Product Datasheet

Bovine Myoglobin ELISA Kit (orb1210600)

Description	Bovine Myoglobin, MB ELISA Kit
Reactivity	Bovine
Range	10-2000ng/mL
Tested Applications	ELISA
Target	MB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,078
Uniprot ID	P02192
Solubility (25°C)	2400ng/mL
Sensitivity	4.47ng/mL
Entrez	280695
Expiration Date	Please enquire.