

Product Datasheet

Rat Ras-related C3 Botulinum Toxin Substrate 1 ELISA Kit (orb1210606)

Description	Rat Ras-related C3 Botulinum Toxin Substrate 1, RAC1 ELISA Kit
Reactivity	Rat
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	RAC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,450
Uniprot ID	Q6RUV5
Solubility (25°C)	96ng/mL
Sensitivity	0.15ng/mL
Entrez	363875
Expiration Date	Please enquire.