

Product Datasheet

Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 6 ELISA Kit (orb1210613)

Description	Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 6, CEACAM6 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	CEACAM6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,195
Uniprot ID	P40199
Solubility (25°C)	64ng/ml
Sensitivity	0.07ng/ml
Entrez	4680
Expiration Date	Please enquire.