

## Product Datasheet

### Mouse Lysosome-associated Membrane Glycoprotein 1 ELISA Kit (orb1210614)

|                            |                                                                    |
|----------------------------|--------------------------------------------------------------------|
| <b>Description</b>         | Mouse Lysosome-associated Membrane Glycoprotein 1, LAMP1 ELISA Kit |
| <b>Reactivity</b>          | Mouse                                                              |
| <b>Range</b>               | 30-9000ng/L                                                        |
| <b>Tested Applications</b> | ELISA                                                              |
| <b>Target</b>              | LAMP1                                                              |
| <b>Storage</b>             | 2-8°C                                                              |
| <b>Note</b>                | For research use only                                              |
| <b>Application notes</b>   | Sandwich                                                           |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                             |
| <b>MW</b>                  | 43,865                                                             |
| <b>Uniprot ID</b>          | <b>P11438</b>                                                      |
| <b>Solubility (25°C)</b>   | 9600ng/L                                                           |
| <b>Sensitivity</b>         | 14.22ng/L                                                          |
| <b>Entrez</b>              | <b>16783</b>                                                       |
| <b>Expiration Date</b>     | Please enquire.                                                    |