

Product Datasheet

Rat Occludin ELISA Kit (orb1210665)

Description	Rat Occludin, OCLN ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	Ocln
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,186
Uniprot ID	Q6P6T5
Solubility (25°C)	3200ng/L
Sensitivity	7.37ng/L
Entrez	83497
Expiration Date	Please enquire.