

Product Datasheet

Mouse Insulin-like Growth Factor-binding Protein Complex Acid Labile Subunit ELISA Kit (orb1210686)

Description	Mouse Insulin-like Growth Factor-binding Protein Complex Acid Labile Subunit, IGFALS ELISA Kit
Reactivity	Mouse
Range	1-400ng/mL
Tested Applications	ELISA
Target	IGFALS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	66,960
Uniprot ID	P70389
Solubility (25°C)	480ng/mL
Sensitivity	0.55ng/mL
Entrez	16005
Expiration Date	Please enquire.