

Product Datasheet

Mouse Inactive Carboxypeptidase-like Protein X2 ELISA Kit (orb1210688)

Description	Mouse Inactive Carboxypeptidase-like Protein X2, CPXM2 ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/mL
Tested Applications	ELISA
Target	CPXM2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	86,963
Uniprot ID	Q9D2L5
Solubility (25°C)	32ng/mL
Sensitivity	0.024ng/mL
Entrez	55987
Expiration Date	Please enquire.