

Product Datasheet

Mouse Icos Ligand ELISA Kit (orb1210691)

Description	Mouse Icos Ligand, ICOSLG ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	ICOSLG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,960
Uniprot ID	Q9JHJ8
Solubility (25°C)	9600ng/L
Sensitivity	14.49ng/L
Entrez	50723
Expiration Date	Please enquire.