

Product Datasheet

Rat Macrophage Colony-stimulating Factor 1 Receptor ELISA Kit (orb1210694)

Description	Rat Macrophage Colony-stimulating Factor 1 Receptor, CSF1R ELISA Kit
Reactivity	Rat
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	CSF1R
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	109,264
Uniprot ID	Q00495
Solubility (25°C)	48ng/ml
Sensitivity	0.044ng/ml
Entrez	307403
Expiration Date	Please enquire.