

Product Datasheet

Mouse Histone Deacetylase Complex Subunit Sap130 ELISA Kit (orb1210697)

Description	Mouse Histone Deacetylase Complex Subunit Sap130, SAP130 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	SAP130
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	111,228
Uniprot ID	Q8BIH0
Solubility (25°C)	6400ng/L
Sensitivity	9.84ng/L
Entrez	269003
Expiration Date	Please enquire.