

Product Datasheet

Mouse Fructose-bisphosphate Aldolase B ELISA Kit (orb1210719)

Description	Mouse Fructose-bisphosphate Aldolase B, ALDOB ELISA Kit
Reactivity	Mouse
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	ALDOB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,507
Uniprot ID	Q91Y97
Solubility (25°C)	96ng/mL
Sensitivity	0.155ng/mL
Entrez	230163
Expiration Date	Please enquire.