

Product Datasheet

Mouse Ezrin ELISA Kit (orb1210724)

Description	Mouse Ezrin, EZR ELISA Kit
Reactivity	Mouse
Range	20-7000ng/L
Tested Applications	ELISA
Target	EZR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,407
Uniprot ID	P26040
Solubility (25°C)	8000ng/L
Sensitivity	9.64ng/L
Entrez	22350
Expiration Date	Please enquire.