

Product Datasheet

Mouse Endothelin-3 ELISA Kit (orb1210727)

Description	Mouse Endothelin-3, EDN3 ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	EDN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,322
Uniprot ID	P48299
Solubility (25°C)	480ng/L
Sensitivity	0.54ng/L
Entrez	13616
Expiration Date	Please enquire.