

Product Datasheet

Mouse Endothelin-2 ELISA Kit (orb1210729)

Description	Mouse Endothelin-2, EDN2 ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	EDN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,538
Uniprot ID	P22389
Solubility (25°C)	480ng/L
Sensitivity	0.56ng/L
Entrez	13615
Expiration Date	Please enquire.