

Product Datasheet

Mouse Endonuclease 8-like 3 ELISA Kit (orb1210732)

Description	Mouse Endonuclease 8-like 3, NEIL3 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	NEIL3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	67,416
Uniprot ID	Q8K203
Solubility (25°C)	4800ng/L
Sensitivity	9.92ng/L
Entrez	234258
Expiration Date	Please enquire.