

Product Datasheet

Mouse Ectonucleotide Pyrophosphatase Phosphodiesterase Family Member 2 ELISA Kit (orb1210735)

Description	Mouse Ectonucleotide Pyrophosphatase Phosphodiesterase Family Member 2, ENPP2 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/mL
Tested Applications	ELISA
Target	ENPP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	98,885
Uniprot ID	Q9R1E6
Solubility (25°C)	4800ng/mL
Sensitivity	9.52ng/mL
Entrez	18606
Expiration Date	Please enquire.