

Product Datasheet

Mouse Disintegrin and Metalloproteinase Domain-containing Protein 17 ELISA Kit (orb1210742)

Description	Mouse Disintegrin and Metalloproteinase Domain-containing Protein 17, ADAM17 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	ADAM17
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	93,056
Uniprot ID	Q9Z0F8
Solubility (25°C)	2400ng/L
Sensitivity	5.23ng/L
Entrez	11491
Expiration Date	Please enquire.