

Product Datasheet

Human Semaphorin 3B ELISA Kit (orb1210748)

Description	Human Semaphorin 3B, SEMA3B ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	SEMA3B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	83,122
Uniprot ID	Q13214
Solubility (25°C)	48ng/ml
Sensitivity	0.046ng/ml
Entrez	7869
Expiration Date	Please enquire.