

Product Datasheet

Mouse Desmosine ELISA Kit (orb1210753)

Description	Mouse Desmosine, DES ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	DES
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,498
Uniprot ID	P31001
Solubility (25°C)	480ng/L
Sensitivity	0.56ng/L
Entrez	13346
Expiration Date	Please enquire.