

Product Datasheet

Mouse C-type Lectin Domain Family 18 Member A ELISA Kit (orb1210768)

Description	Mouse C-type Lectin Domain Family 18 Member A, CLEC18A ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	CLEC18A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,585
Uniprot ID	Q7TSQ1
Solubility (25°C)	48ng/mL
Sensitivity	0.056ng/mL
Entrez	353287
Expiration Date	Please enquire.