

Product Datasheet

Mouse Complement C1Q Subcomponent Subunit A ELISA Kit (orb1210776)

Description	Mouse Complement C1Q Subcomponent Subunit A, C1QA ELISA Kit
Reactivity	Mouse
Range	0.05-20ug/mL
Tested Applications	ELISA
Target	C1QA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,974
Uniprot ID	P98086
Solubility (25°C)	24µg/mL
Sensitivity	0.023µg/mL
Entrez	12259
Expiration Date	Please enquire.