

Product Datasheet

Mouse Cold-inducible RNA-binding Protein ELISA Kit (orb1210780)

Description	Mouse Cold-inducible RNA-binding Protein, CIRBP ELISA Kit
Reactivity	Mouse
Range	5-1000ng/L
Tested Applications	ELISA
Target	CIRBP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,607
Uniprot ID	P60824
Solubility (25°C)	1280ng/L
Sensitivity	2.24ng/L
Entrez	12696
Expiration Date	Please enquire.