

## Product Datasheet

### Rat Creatine Kinase M-type ELISA Kit (orb1210790)

|                                |                                              |
|--------------------------------|----------------------------------------------|
| <b>Description</b>             | Rat Creatine Kinase M-type,<br>CKM ELISA Kit |
| <b>Reactivity</b>              | Rat                                          |
| <b>Range</b>                   | 0.6-180mU/mL                                 |
| <b>Tested<br/>Applications</b> | ELISA                                        |
| <b>Target</b>                  | CKM                                          |
| <b>Storage</b>                 | 2-8°C                                        |
| <b>Note</b>                    | For research use only                        |
| <b>Application notes</b>       | Sandwich                                     |
| <b>Sample Types</b>            | Serum, plasma, cell culture<br>supernates    |
| <b>MW</b>                      | 43,045                                       |
| <b>Uniprot ID</b>              | <b>P00564</b>                                |
| <b>Solubility (25°C)</b>       | 192mU/mL                                     |
| <b>Sensitivity</b>             | 0.32mU/mL                                    |
| <b>Entrez</b>                  | <b>24265</b>                                 |
| <b>Expiration Date</b>         | Please enquire.                              |