

Product Datasheet

Mouse Cd9 Antigen ELISA Kit (orb1210795)

Description	Mouse Cd9 Antigen, CD9 ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	CD9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,258
Uniprot ID	P40240
Solubility (25°C)	640ng/L
Sensitivity	1.27ng/L
Entrez	12527
Expiration Date	Please enquire.