

Product Datasheet

Mouse Carnitine O-palmitoyltransferase 1 ELISA Kit (orb1210804)

Description	Mouse Carnitine O-palmitoyltransferase 1, CPT1A ELISA Kit
Reactivity	Mouse
Range	2-700ng/L
Tested Applications	ELISA
Target	CPT1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	88,251
Uniprot ID	P97742
Solubility (25°C)	800ng/L
Sensitivity	1.16ng/L
Entrez	12894
Expiration Date	Please enquire.