

Product Datasheet

Mouse Carbonic Anhydrase 3 ELISA Kit (orb1210812)

Description	Mouse Carbonic Anhydrase 3, CA3 ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	CA3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,366
Uniprot ID	P16015
Solubility (25°C)	640ng/L
Sensitivity	1.27ng/L
Entrez	12350
Expiration Date	Please enquire.