

Product Datasheet

Rat C-C Motif Chemokine 6 ELISA Kit (orb1210815)

Description	Rat C-C Motif Chemokine 6, CCL6 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CCL6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,652
Uniprot ID	Q68FP3
Solubility (25°C)	24ng/ml
Sensitivity	0.025ng/ml
Entrez	287910
Expiration Date	Please enquire.