


Product Datasheet

Rat Caspase-1 ELISA Kit (orb1210819)

Description	Rat Caspase-1, CASP1 ELISA Kit
Reactivity	Rat
Range	3-900ng/L
Tested Applications	ELISA
Target	CASP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,576
Uniprot ID	P43527
Solubility (25°C)	960ng/L
Sensitivity	1.53ng/L
Entrez	25166
Expiration Date	Please enquire.



Standard Curve of Rat Caspase-1 ELISA Ki...