

Product Datasheet

Rat Calcium Calmodulin-dependent Protein Kinase Type II Subunit Delta ELISA Kit (orb1210824)

Description	Rat Calcium Calmodulin-dependent Protein Kinase Type II Subunit Delta, CAMK2D ELISA Kit
Reactivity	Rat
Range	20-4500ng/L
Tested Applications	ELISA
Target	CAMK2D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	60,081
Uniprot ID	P15791
Solubility (25°C)	4800ng/L
Sensitivity	9.12ng/L
Entrez	24246
Expiration Date	Please enquire.