

Product Datasheet

Rat Beta-1 Adrenergic Receptor ELISA Kit (orb1210830)

Description	Rat Beta-1 Adrenergic Receptor, ADRB1 ELISA Kit
Reactivity	Rat
Range	20-4500ng/L
Tested Applications	ELISA
Target	ADRB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,472
Uniprot ID	P18090
Solubility (25°C)	4800ng/L
Sensitivity	9.05ng/L
Entrez	24925
Expiration Date	Please enquire.